

Oxybaby® Control Center



OBCC Software


Software for gas analyser OXYBABY® / PA

Windows-Software for OXYBABY® 6.0, OXYBABY® P 6.0, OXYBABY® Med and PA 7.0. OBCC is the ideal complement for your Gasanalyser - with modern surface, clearly structured and conveniently operable.

The application offers rapid graphical analysis of measured values. Simplify a system comprehensive data management and offers an uncomplicated export function for the complete digital storage.

Benefits

- Modern, clearly structured user interface
- Intuitive operation
- Simple set up and administration of master data (products, user, packaging lines)
- Quality control by means of data recording and storing with date and time
- The graphical presentation of measured values offers a fast overview
- Knowledge by data analysis tool with filter options
- Completely digital documentation by simple CSV export function
- 1 click generate PDF report for assurance of proof (in addition to the CSV export)
- Time-saving by convenient data management and easy transfer of master data to new devices
- Available menu languages: English, German, French, Italian, Spanish, Portuguese, Polish, Russian, Finnish, Swedish, Turkish, Romanian, Hungarian



System requirements

Operating system	WINDOWS® 7 or above
Hardware	1 free USB interface, graphic resolution min. 1024 x 768 Pixel

HTK Hamburg GmbH
 Oehleckerring 32 • 22419 Hamburg
 Phone: +49 (0)40 - 600 38 38 - 0 • Fax: +49 (0)40 - 600 38 38 - 99
info@htk-hamburg.com



© Copyright 2019 - All contents of this document, in particular Texts, photographs and graphics are protected by copyright. All rights, including reproduction, publication, processing and translation are reserved, HTK Hamburg GmbH. Please contact HTK Hamburg GmbH if you would like to use the contents of this document.

Rev.III_022021_Oxybaby Control center_deu • Subject to change



We control GASES - since 1978

